



LuxCreo LuxAero

LuxAero is LuxCreo's dedicated odor-control solution for LuxOven Pro, designed to capture and filter exhaust directly from the source and significantly reduce noticeable odors during post-processing.

LuxAero delivers a cleaner, more comfortable post-processing experience for dental labs and clinics.

LuxAero



SYSTEM SPECIFICATIONS

System Airflow	165m ³ /h
Exterior Dimension Uncrated	285 x 195 x 280 mm (11.2 x 7.7 x 11.0 in)
Exterior Dimension Crated	360 x 320 x 390 mm (14.2 x 12.6 x 15.4 in)
Uncrated Weight	6.4 kg (14.1 lbs)
Crated Weight	7.7 kg (16.98 lbs)
Filtration Efficiency	0.3µm / 99.97%
Noise Level	<50dB

SUGGESTED OPERATING ENVIRONMENT

Temperature	22-32 °C (72-90 °F)
Humidity (RH)	<80%
Power Requirement	AC 100~240V, 50/60 Hz, 0.90~0.38A